1) PRODUCT DESCRIPTION

LIQUEEL GD is a sterile, absorbent, self-adhesive, soft wound dressing made of polyurethane foam with a silicone coating. It consists of a soft wound-contact layer, a highly absorbent polyurethane foam layer and a vapour-permeable outer layer made of polyurethane film.

2) PRODUCT FEATURES

In the presence of exudate, LIQUEEL GD helps to maintain a moist wound environment in order to support the natural healing process.

3) INDICATIONS

LIQUEEL GD is indicated for the treatment of lightly to moderately exuding wounds, such as

- pressure sores
- venous and arterial leg ulcers
- non-infected diabetic foot ulcers
- first- and second-degree scalds/burns
- abrasions, lacerations, and small cuts

4) CONTRAINDICATIONS

LIQUEEL GD is contraindicated for infected wounds, deep wounds, heavily bleeding wounds, and greater than second-degree burns. LIQUEEL GD should not be used on patients with a known hypersensitivity to polyurethane.

5) SAFETY INFORMATION

The wound dressing is sterile as long as the packaging is not damaged or opened prior to use.

The product may only be used under these conditions.

LIQUEEL GD is only intended for single use and must not be used after its expiry date.

5.1) WARNINGS

- Do not use if pouch is damaged or opened.
- Do not re-use.

Re-use of single-use devices creates a potential risk to the patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the product may lead to injury, illness or death of the patient.

LIQUEEL GD must not be used in combination with oxidising agents such as hydrogen peroxide.

6) INFORMATION ON USE

LIQUEEL GD is easy to apply and does not require the use of sterile instruments. The interval between dressing changes depends on the individual wound condition.

In the event of intolerance to the wound dressing, it should be removed and replaced with a tolerated alternative.

6.1) PREPARATION

- Visually inspect pouch for damage prior to opening.
- Clean the wound with a sterile wound irrigation solution such as saline or Ringer's solution, and sterile swabs and/or pads.
- Carefully dry the skin surrounding the wound.

6.2) APPLYING THE DRESSING

- Remove the sterile wound dressing from the pouch using local guidelines & best practices for handling sterile wound dressings.
- Remove the smaller film section from the back of the wound dressing.
- Align the middle of the wound dressing with the centre of the wound, and apply the bottom part of the wound dressing to the intact skin.
- If the shape or size of the wound dressing needs to be modified, tear along the perforations, as required. If no modification is required, remove the remaining film from the back and continue with step 7.
- Provisionally modify the wound dressing to suit the size of the wound and the anatomy by detaching one or more perforated sections.
- After the provisional modification, remove the remaining film from the side of the dressing that faces the wound.

- After removing the protective film, centre the dressing over the wound and apply.

An appropriate secondary dressing may be used, e.g. a PU film.

6.3) CHANGING THE DRESSING

You should replace dressing LIQUEEL GD if it becomes soiled, saturated or if exudate/ drainage is observed or adhesion is compromised. Otherwise, replace dressing LIQUEEL GD per established facility protocol.

- Remove LIQUEEL GD gently
- Follow steps 6.1 - a) to 6.2 - 7 to apply a new dressing.

7) PRECAUTIONS AND SPECIAL INSTRUCTIONS

LIQUEEL GD can be used under compression therapy. Because of its good adhesive properties, LIQUEEL GD should be removed with care if the skin is sensitive (e.g. parchment skin).

As part of the normal healing process, the wound will shed dead tissue. This is desirable and beneficial to wound healing. As a result, the wound may initially appear to be getting larger.

The frequency of dressing changes will depend on the quantity of exudate and the localised wound healing factors, and should be adjusted accordingly.

8) STORAGE/STERILISATION/SYMBOLS ON THE LABEL

Store the dressings away from direct sunlight, in a dry place and at room temperature.

9) DISPOSAL

The silicone foam dressings should be disposed of based on the hospital or healthcare professional advice.

DO NOT USE THE PRODUCT IF THE POUCH IS DAMAGED OR OPEN.

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 Sterile product; sterilised with ethylene oxide

Caution

Follow the Instructions for Use.

Single sterile barrier system.

Do not re-sterilise!

Single use only.
Do not reuse.

Product is sterile as long as sterile packaging is not open or damaged.

Keep dry.


Keep away from sunlight.

Avery Dennison Medical Ltd.
 IDA Business Park, Ballinalee Road
 Longford, N39 DX73. Ireland


REV_NO_LQGD_20201129_1225_EN


1)

Remove the sterile wound dressing from the pouch using local guidelines & best practices for handling sterile wound dressings.


2)

Remove the smaller film section from the back of the wound dressing.


3)

Align the middle of the wound dressing with the centre of the wound, and apply the bottom part of the wound dressing to the intact skin.


4)

If the shape or size of the wound dressing needs to be modified, tear along the perforations, as required. If no modification is required, remove the remaining film from the back and continue with step 7.


5)

Provisionally modify the wound dressing to suit the size of the wound and the anatomy by detaching one or more perforated sections.

6)

After the provisional modification, remove the remaining film from the side of the dressing that faces the wound.

7)

After removing the protective film, centre the dressing over the wound and apply.

REV_NO_LQGD_20201129_1225_EN

		PREPRESS NOTES			Proof Revision: 05	Comments:
<input checked="" type="checkbox"/> OK to Proceed	Item Number:	100045A - LIQUEEL GD	PDF Page:	1 of 2	Operator:	Doreen Mulroy
<input type="checkbox"/> Make Corrections and Proceed	Description:	N/A	Dimensions:	420 x 297mm folded to 105 x 148mm	Resolution Proof:	300dpi
<input type="checkbox"/> Make Corrections and Submit Another Proof	Client Name:	Avery Dennison Medical	Inks:	CMYK	Date:	1/12/2020
Signed  Date 03.12.2020			Page Count:	2	Time:	13:03

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