FRONT (final fold)

BACK (final fold)

90 mm

90 mm

90 mm

45 mm

90 mm

270 mm

BD ChloraShield™ IV

BD ChloraShield™ IV Surround

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REF 410111, REF 410122

14064A
THT-PKG-000185-A
5-775111-C

ADVERSE REACTIONS SUCH AS IRRITATIONS, SENSITIZATION AND GENERALIZED ALLERGIC REACTIONS HAVE BEEN REPORTED WITH THE USE OF CHLORHEXIDINE GLUCONATE.

- STOP USING THE BD CHLORASHIELD™ IV DRESSING WITH CHG ANTIMICROBIAL IF THE PATIENT EXPERIENCES SYMPTOMS SUCH AS: WHEEZING OR DIFFICULTY BREATHING, SWELLING FACE, HIVES THAT CAN QUICKLY PROGRESS TO MORE SERIOUS SYMPTOMS, SEVERE RASH OR SHOCK (WHICH IS A LIFE THREATENING CONDITION THAT OCCURS WHEN THE BODY IS NOT GETTING ENOUGH BLOOD FLOW). THESE SYMPTOMS ARE ASSOCIATED WITH AN ALLERGIC REACTION TO CHG.
- Literature searches found hypersensitivity reactions associated with the topical use of chlorhexidine gluconate have been reported in several countries. Caution should be used when using chlorhexidine containing preparations, and the patient should be observed for the possibility of hypersensitivity reactions.
- IF ALLERGIC REACTIONS OCCUR, DISCONTINUE USE IMMEDIATELY AND, IF SEVERE, CONTACT A PHYSICIAN.
- INTENDED FOR SINGLE PATIENT USE ONLY.
- DO NOT REUSE. AS WITH ALL ADHESIVE-BASED PRODUCTS, ADHESIVE EFFECTIVENESS AND FUNCTIONALITY CAN DECLINE AFTER THE FIRST USE AND THE PRODUCT WILL NOT PERFORM AS SPECIFIED. REUSE MAY LEAD TO INFECTION OR OTHER ILLNESS/INJURY.
- STORE IN COOL, DRY PLACE. KEEP AWAY FROM DIRECT SUNLIGHT. DO NOT STORE ABOVE 35°C/95°F. DO NOT USE IF PACKAGING IS DAMAGED. DO NOT RESTERILIZE.

Instructions for Use:
Failure to follow the manufacturer's instructions may result in complications including skin irritation and/or maceration. If any damage to the packaging is observed, the sterility of the dressing may be compromised. Do not use the dressing and contact the BD Customer Support team at 1.800.323.9088.

Site Preparation:

- Prepare the site in accordance to facility protocol.
- Stabilize any active bleeding before applying the dressing.
- Ensure the skin is clean, dry and free from detergent, lotion and soap residues.
- Ensure all skin preparations and protectants dry completely before applying the dressing. This will promote optimal adhesion to the skin and help minimize any skin irritation.
- Not all skin types or conditions will allow for full adhesion, therefore, adhesive performance may vary between patients.

Dressing Application:
75mm X 60mm (REF 410111)

- Remove the dressing from the sterile package using aseptic technique.
- Starting at the orange end, peel the bottom liner (marked 1) from the back of the dressing, fully exposing the adhesive surface.
- Hold the dressing by the extended adhesive-free orange tab and center the dressing over the catheter site. Provide at least a one inch border around the catheter insertion site.
- Apply the adhesive side onto the skin and smooth dressing into place.
- Apply pressure to the edge of the dressing all the way around the perimeter to minimize potential for edge lift.
- Remove the top cover film by lifting up on the white tab (marked 2) and slowly pull in the direction of the arrows to fully remove the white tab and cover film. Rub down any edges of the dressing that may not be fully adhered.
- Break off the extended orange tab at the dotted line and discard it.
- Document the dressing change in accordance with facility protocols. The data strip on the dressing may be utilized for this purpose.

75mm X 60mm with nonwoven border and slit (REF 410122)

- Remove the dressing from the sterile package using aseptic technique.
- Starting at the orange end, peel the bottom liner (marked 1) from the back of the dressing until the liner reaches a fold, exposing part of the adhesive surface.
- Hold the dressing by the extended orange tab with one hand and, using the other hand, hold the partially removed liner back along with the other half of the dressing to keep the exposed adhesive area clear for visibility during placement.
- Apply the exposed adhesive area over the catheter insertion site and smooth dressing into place.
- Remove the remainder of the bottom liner and smooth the rest of the dressing around and under the catheter. Ensure the edges of the slit are joined together under the catheter to fully surround the insertion site.
- Break off the extended orange tab at the dotted line and discard it.
- Document the dressing change in accordance with facility protocols. The data strip on the dressing may be utilized for this purpose.

Site Care:

- The site should be observed at least daily for signs of infection or other complications. If infection is suspected, remove the dressing, inspect the site directly, and determine the appropriate medical intervention. Infection may be signaled by fever, pain, redness, swelling, or unusual discharge or odor.
- Change the dressing as necessary, in accordance with facility protocols. At a minimum, dressings should be changed at least every 7 days, per current CDC recommendations. Total dressings use should not exceed 30 consecutive days, regardless of how often the individual dressing is changed. Change the dressing more frequently if the site has high levels of exudate or if the dressing integrity has become compromised.

3. The dressing should be changed sooner than 7 days if:

- the dressing becomes loose, soiled or compromised in any way
- the site is obscured or no longer visible
- there is visible drainage underneath the dressing

Removal:

- Gently grasp the corner of the dressing and slowly peel it from the skin in the direction of hair growth.
- To minimize chance of skin trauma, remove dressing using the low and slow technique, pulling the dressing in the direction of the catheter site.
- When it is desired to maintain catheter in place during a dressing change, direct pressure should be maintained through the dressing onto the catheter hub while the dressing is loosened around and over part of the catheter hub.
- Then, once most of the dressing has been removed using the low and slow technique, apply direct pressure to the exposed portion of the catheter hub as the dressing is removed from the remainder of the hub.

If needed, medical adhesive solvents can be used to help remove the dressing.

4. Discard in accordance with your facility protocol.

Additional Information:

Product REF	410111
Product Name	BD ChloraShield™ IV Dressing with Chlorhexidine Gluconate (CHG) Antimicrobial
Dressing Size	75 mm x 60 mm (3.0 in x 2.4 in)
Quantity per Box	100

Product REF	410122
Product Name	BD ChloraShield™ IV Surround Dressing with Chlorhexidine Gluconate (CHG) Antimicrobial
Dressing Size	75 mm x 60 mm (3.0 in x 2.4 in)
Quantity per Box	100

For questions or comments, contact the BD Customer Support at 1.800.323.9088.

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An Avery Dennison Business.

THE INFORMATION LISTED BELOW MUST APPEAR ON ALL VENDOR PROOFS (COLORKEYS, MATCHPRINTS/CHROMALINS) OR THEY WILL BE REJECTED							
AW	DESCRIPTION	SPEC. NO.	PAGE	PMS COLORS	PROCESS MATCH COLORS	SPECIAL INSTRUCTIONS	APPROVAL SIGNATURE / DATE
THT-PKG-000185-A, 14064A	BD ChloraShield IFU, IV Gen 1 and Gen 2 U.S. Region		1 of 2	■ PANTONE Black C	NONE	Red die lines do not print. Dimensions do not print.	

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270 mm

An Avery Dennison Business.

THE INFORMATION LISTED BELOW MUST APPEAR ON ALL VENDOR PROOFS (COLORKEYS, MATCHPRINTS/CHROMALINS) OR THEY WILL BE REJECTED							
AW	DESCRIPTION	SPEC. NO.	PAGE	PMS COLORS	PROCESS MATCH COLORS	SPECIAL INSTRUCTIONS	APPROVAL SIGNATURE / DATE
THT-PKG-000185-A, 14064A	BD ChloraShield IFU, IV Gen 1 and Gen 2 U.S. Region		2 of 2	■ PANTONE Black C	NONE	Red die lines do not print. Dimensions do not print.	

DCN IP80158 5-775111-C